# INSTRUCTION MANUAL

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## RELAXMUS

### Fascia Mini Massager

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### CONTAINS:

1. RELAXMUS Fascia mini massager
2. Massager heads x 4
3. Carrying case
4. Hanging strap
5. USB-C to USB-A Cable (50cm)
6. Instruction manual and warranty card

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In the charging state, when the power indicator light is green, it means that it is charging. The green indicator light is off when fully charged. When charging, please make sure that the massager is turned off. 充電狀態下，當電量指示燈為綠色時，表示正在充電，充滿電綠色指示燈熄滅，充電時請確保按摩器為關機狀態。

The runtime of the device is about 4 hours (max.), depending on the vibration frequency and applied pressure. 設備充電後運行時間約4個小時(最多)，時間取決於使用期間的振動頻率和施加的壓力。

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Quick Installation 快速安裝 / 快速安装

## INSTRUCTION FOR USE

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In the power off state, press and hold the ① button on the panel for 2 seconds to turn on and enter level 1, short press the ① button on the panel, the level will be switched between level 1 to level 5; Press and hold the ① button for 2 seconds when the device is in the level 1-5, it will enter standby mode; When the product is in standby mode, press and hold ① for 2 seconds to turn off the power.

After the product operates for 10 minutes, it will automatically enter standby mode; LED lights indicate the corresponding vibration frequency level, 1-level is the lowest, 5-level is the highest; When the device is in standby mode, it will automatically power off after 20 seconds idle time.

Apply the desired pressure at the preset vibration frequency to massage the corresponding body parts.

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About the recommended groups of people and parts of use:

- Only for healthy adults except the elderly;
- Use only on the soft tissues of the body when not causing pain or discomfort;
- Do not use the product on the head or any hard or bony part of the body;
- Bruising may occur with either moderate or applied pressure, please check the massaged area frequently and stop immediately at the first sign of pain or discomfort;
- To be used only on dry, cleansed skin surfaces on body parts, gently Press and move on the skin;
- Keep fingers, hair or other body parts away from the shaft and back of the massage head, as the possibility of being caught;
- Do not use near bony areas with less muscle coverage such as: tibia, back of foot, back of hand;
- Do not use on and around sensitive areas: eyes, teeth, breast augmentation, reproductive organs, thetics, etc.;
- Do not use the product on any possible fracture site;

- If pain persists, discontinue use immediately and consult a physician.
- This product is not a medical device and is not recommended for medical use and is not intended to be used in hospitals for medical use;
- Please check the main unit and connecting wire before use, if there is any damage, stop using it to avoid danger.

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Copyright © 2024 MOMAX Technology Ltd. All rights reserved. MOMAX is registered trademark of MOMAX Technology Ltd. All other trademarks are the property of their respective owners. MOMAX是MOMAX Technology Ltd. 保留所有權利。所有其他商標均為其各自擁有者的財產。 版權所有 ©2024 MOMAX Technology Ltd. 保留所有權利。 MOMAX是MOMAX Technology Ltd. 的註冊商標。所有其他商標均為其各自擁有者的財產。